

## LASOS<sup>®</sup> System series



Multi Color Systems

**MCS 4 series, MCSL 4 series and MCS 6 series**

Multi-wavelength laser modules and beam combiners

- Different designs for compact and flexible solutions
- Benchtop devices and fully installable sub-systems
- Free beam or fiber output
- Different options for beam collimation and combining



## MCS 4 series / MCSL 4 series



Combines up to four wavelengths into one or two fiber outputs

RGB-broad band fiber with optional separate IR filter

Customized wavelength and output power combinations

Standard fiber connectors or collimated beam output



### ADDITIONAL CHARACTERISTICS

Each wavelength can separately be modulated

Passive cooling via device plate

Interface for Mini-USB, Input/Output and 12 V power supply





### AVAILABLE WAVELENGTHS [nm]

Please refer to LASOS DPSS and LDM laser series for available wavelengths

Possible combination within the system: 1 DPSS and 3 LDM or 4 LDM

### DIMENSIONS [mm]



## MCS 6 series



Delivers up to six wavelengths with fiber output from a standard 19" housing

Customized wavelength and output power combinations

Standard fiber connectors or collimated beam output


Optional beam combining unit with separately switchable laser lines



### ADDITIONAL CHARACTERISTICS

Each wavelength is individually changeable und upgradeable


Laser lines are separately fan-cooled controlled by the ambient temperature



### AVAILABLE WAVELENGTHS [nm]

Please refer to LASOS DPSS and LDM laser series for available wavelengths

### DIMENSIONS [mm]



### FEATURES

- Up to 100 mW output depending on the wavelength
- Customized configurations and modifications available
- Standard fiber connectors or collimated beam output
- Fixed or separable fiber coupling
- Single mode, multimode or polarization maintaining fibers

### TECHNICAL DATA

- Operational temperature 15°C - 35°C
- Output power stability  $\pm 2\%$



LASOS Lasertechnik GmbH  
Franz-Loewen-Strasse 2  
07745 Jena  
Germany

Phone: +49 3641 2944 0  
Fax: +49 3641 2944 300  
info@lasos.com

[www.lasos.com](http://www.lasos.com)